Zeiss AxioObserver with ApoTome Zen software

User Guide

LSU Health-Shreveport

Research Core Facility (RCF)

Chaowei Shang

July 30th 2021



If the room lights are not on when you arrive to use the Axioobserver, **DO NOT TURN THEM ON**. Please be courteous of those using the other microscopes, since they may have light-sensitive slides, dishes, or plates on the microscope stages, which could be damaged by the room lights. Use a flashlight if needed.

Check whether the stage **holder** is properly in place. If not, please contact the Research Specialist. Make sure the **table is clean**. No food or water is allowed in the microscope room.

Before you start, please check the objective turret has the 10x objective in place. This is a safety precaution. If the 10x objective is not in place, please clean the current objective, drop the objective to the bottom, and switch it to the 10x objective. Please also inform the specialist that the objective is not positioned to 10x before you start.



Stage holder

START UP THE SYSTEM

Please sign in on the login sheet. The format for the date is MM/DD/YYYY, eg. 05/25/2015. Write legibly for this is how the staff knows who has been using the microscope. Check in and Check out time is written in 24 hour format, e.g. 1pm is 13. The format for writing the bulb time is ###h ## m, eg. 100h 25m. (You will not be able to write the bulb time until you turn on the power supply.)

The computer should be on. Log in with your **LSUHSC ID and password**. Make sure the domain is set to **LSUMC-MASTER**, selected from the drop-down menu. If you only wish to access documents/files, then login and do not turn on the mercury lamp power supply or the other buttons. Please indicate this on the sign in sheet.

Please sign in on the login sheet. The format for the date is MM/DD/YYYY, eg. 05/25/2015. Write legibly for this is how the staff knows who has been using the microscope. Check in and Check out time is written in 24 hour format, e.g. 1pm is 13. The format for writing the bulb time is ###h ## m, eg. 100h 25m. (You will not be able to write the bulb time until you turn on the power supply.)

Our system has the startup labeled numerically with yellow stickers.

- 1. Turn on the mercury lamp power source. The green power light and the green lamp light will come on. (This is where you record the bulb time.) If the lamp light indicator does not light up, please contact Core facility staff. Once turned on, the lamp should remain on for a minimum of 30 minutes. Do not turn the lamp off if another user is scheduled within 2 hours after your session.
- 2. Turn on the Power Supply 2 by flipping the switch which will have a green light turn on.
- 3. Turn on the microscope (number 3). Button located at the left side of the microscope. This will turn on the Touch Screen Controller, and please wait for this display to fully load before proceeding.
- 4. Turn on the Apotome power supply labeled with a yellow 4 sticker if you are going to use apotome.



Touch Screen Controller









TOUCH SCREEN CONTROLLER

Focus knob: Find the focus (move objective up and down) Two "wheels" for coarse and fine adjustments and controls the Z-direction of movement. **Joystick**: Find the sample (move the stage). The joystick controls the X and Y direction of the stage

An icon is selected when it is in white and not selected when in blue.

Home screen information:

- 1. **Objective**: Description of the current objective in the light path along with N.A.
- 2. Resolution: How fine of detail the current image will appear
- 3. Total Magnification: magnification of objective * eyepiece magnification (10x)
- 4.Reflector: Displays which position of the reflector turret is in the light path and what type of cube is in that position.
- 5.VIS: Displays how much of the light from the sample is being sent to the eyepiece
- 6.Sideport: Displays how much of the light from the sample is sent to the camera/computer
- 7. Z-position: The bigger the number the higher the objective, the smaller the number the lower the objective.

8: X, Y-position: sample position coordinates.

- 9: Load Position (on the right side of the screen): Press this will drop the objective to the bottom.
- 10: TL Illumination: Transmitted light on or off
- 11: RL Illumination: Fluorescent light on or off











On the left of the screen, push the icon labeled Microscope. A new window will pull up in which you can select Objectives, Reflector, and Light path.



In the Microscope screen, press the **Objectives** tab to enter the objectives screen. Also, ensure that you are on the **Control** display by checking the **Control** icon on the left sidebar. The **Objectives** displayed tell you which one is an air objective versus an oil objective. If you need to switch between different objectives that use different immersion, the software and TFT display will inform you that you should be prepared for a different immersion objective

Press Reflector on the top of the microscope screen and you will enter the Reflector screen.

DIC: Transmitted light

34 BFP: DAPI channel (EX 390/20, EM 460/60) **38 HE GFP: Green GFP channel** (EX 470/40, EM 525/50) **31 AF 568: Red/ DsRed channel** (EX 565/30, EM 620/60) **17 AF 488: FITC channel** (EX 485/20, EM 540/25) **50 Cy5: 647 far red channel** (EX 640/30, EM 690/50)



Reflector Screen

Press Light Path on the top of the microscope screen and you will enter the Light Path screen



Sideport L: Directs light to go through the camera and images will show up on the computer.

Eyepiece/VIS: Directs light to go through Eyepiece/Vision, so you can observer through eyepiece.

Press the square next to the Camera and you will have 3 options (100% camera, 100% eye, and 50% Camera & 50% eye). Do not press the square next to Eye, it won't change anything

Everything on the touch screen can be controlled by the software. It is preferred to control everything through the software, but below shows you briefly how to set up for imaging if you do not use the software.

Set up for eyepiece observation

- 1. Select an objective through the Objective screen.
- 2. Go to the Reflector screen and select a channel. Then turn off TL, and turn on RL
- 3. Go to the Light path screen and select 100% light goes to Eyes.
- 4. Locate your sample and focus.
- 5. When done, press RL illumination off to turn off the light (Very important!!)





Set up for the software to acquire an image

After turning off RL illumination, switch light path to 100% to sideport



Start up the Zen software

To acquire images, all hardware MUST be turned on before you start up the software.

1. On the desktop, click the Zen 3.3 (blue edition) icon to launch the software.

2. Then select Zen pro for image acquisition. If you only want to look at an image, or do image analysis, click ZEN image processing.

3. Next, click Skip Calibration when you see this window pop up.





If you only want to look at an image, or do image analysis using ZEN image processing. You do not need to turn on the microscope hardware.

File Edit View Acquisition Graphi	ics Tools Window Help
Scaling	5 A
Locate Acquisition Processing	Analysis Applications
Transmitted Light Off On Re	flected Light Off On
Favorites Configure	
DAPI GFP FITC	DsRed DIC
Eyes Camera	
Find Focus Set Exposure Live	Continuous Snap
Active Camera AxioCam MR R3	
Microscope Control	Show All
Camera	
Exposure Time	
Time (50.0 🗘 ms 🔻
Auto Exposure Intensity 10	0 % 🗘 Set Exposure
Spot Meter / Focus ROI	
Binning 1x1 - Resolution	on 1388x1040 -
Binning-independent Brightness	
Acquisition ROI	
• Gain	
Post Processing	
 Settings 	
Mode	
Trigger Control	
Model Specific	
Software Autofocus	Show All
E Movie Recorder	
🕨 📖 Live Panorama	

Locate section



When software turns on, by default, it displays the Locate tab/section (on the left)

The **Locate** section allows you to turn on the fluorescent light and view the sample through eyepieces.

The Objectives section on the right allows you to change objectives.

Set up eyepiece observation through software

1. In the Locate section, click to select a fluorescent/Bright field channel , (DAPI, GFP, FITC, DsRed, DIC). This will turn on the light on the microscope.

2. When light is turned on, the lightpath automatically switch to eyes, you do not need to press anything to make this happen.

3. Then use the joystick to locate your sample. Use the focus knob to focus on your sample.

4. When done, turn off the fluorescent light either through pressing the RL/TL Illumination Off on the touch screen controller, or turn off the light through pressing the Off button on the software.

It is VERY IMPORTANT to turn off the light after observation, so light does not photobleach your sample.

Scaling Cocate Acquisition Processing Analysis Applications Transmitted Light Off On Reflected Light Off Off On Reflected Light Off Off On Reflected Light Off On On Reflected Light Off On On On Reflected Light On <po< th=""></po<>
Image: Control Arive Camera Active Cam
Transmitted Light Off On Reflected Light Off On Favorites Configure OAN GFP FIC Defked DIC Eyes Camera AF Image: Continuous Snap Active Camera AxioCam MR R3 Microscope Control Show All Image: Comera Show All Exposure Time 0.0 ms Time 60.0 ms Auto Exposure Intensity 10% Spot Meter / Focus ROI Binning 1X1
Favorites Configure DAP GFP FITC DRRedi DIC Eyes Camera Image: Continuous Snap AF Image: Continuous Snap Image: Continuous Snap Active Camera AxioCam MR R3 Image: Control Show All Image: Control Image: Microscope Control Show All Image: Control Show All Image: Control Image: Comera Go.0 Image: Control Show All Image: Control Show All Image: Control Show Line Image: Control Show All Image: Control Image: Control Show All Image: Control
DARI GFP FITC DaRed DIC Eyes Camera AF Image: Continuous Snap Active Camera AxioCam MR R3 Active Camera Show All Image: Camera Show All Image: Camera Show All Exposure Time 60.0 Time 60.0 Spot Meter / Focus ROI Binning IXT Resolution 1388x1040
Eyes Camera AF Image: Common set of the
AF Find Focus Set Exposure Live Continuous Snap Active Camera AxioCam MR R3 Microscope Control Show All Microscope Control Show All Show All Exposure Time Time 60.0 rms Auto Exposure Intensity 100 % Set Exposure Spot Meter / Focus ROI Binning 1x1 Resolution 1388x1040
Find Focus Set Exposure Live Continuous Snap Active Camera Active Camera Show All Image: Show All Image: Show All Exposure Time Time 60.0 Image: Show All Image: Auto Exposure Intensity 100 % Set Exposure Spot Meter / Focus ROI Binning Ix1 Resolution 1388x1040
Active Camera AxioCam MR R3
 Microscope Control Show All Camera Show All Show All Exposure Time Go.0 ms Go.0 ms Auto Exposure Intensity 100 % Set Exposure Spot Meter / Focus ROI Binning Ix1 Resolution I388x1040
Exposure Time Time 60.0 Auto Exposure Intensity Spot Meter / Focus ROI Binning 1x1
Exposure Time Time 60.0 ms Auto Exposure Intensity 100 % Set Exposure Spot Meter / Focus ROI Binning 1x1 Resolution 1388x1040 v
Time 60.0 ms Auto Exposure Intensity 100 % Set Exposure Spot Meter / Focus ROI Binning 1x1 Resolution 1388x1040
Auto Exposure Intensity 100 % Set Exposure Spot Meter / Focus ROI Binning 1x1 Resolution 1388x1040
Spot Meter / Focus ROI Binning 1x1 Resolution 1388x1040
Binning 1x1 Resolution 1388x1040
Binning-independent Brightness
Acquisition ROI
O Gain
Post Processing
O Settings

View and acquire image through the software:

- Click on Acquisition to switch to the Acquisition section.
- The middle window is where the live image will display.
- The bottom is the histogram and image tools.



The **Acquisition** section allows you to view and capture images through the software.

- 1. Click on the drop down menu as shown in the picture. Select All Channels. This will bring up a list of all the fluorescent channels.
- Check the box in front of the channel to select the channels you need.

		Locate
own in the I bring up a l	ist	Experim User Do All Chai Apotor
Acquisitio	on Mode	Show All
	Acquisition	Show All
		gr
🗸 DAPI		Ref. 🗖 🗖 ion
EGFP		
🗆 FITC		
DsRed		- -
🗆 Су5		
	Focus Ref.	* - * -
Lightsource		
DAPI		
Camera	AxioCam MR R3	
Exposure	Auto Exposure	Set Exposure
Time		- <u>300.0</u> ; ms ▼



To start a live view of the image

- 1. Click on the selected channel so that it turns light grey color.
- 2. Then click on Live, and you will see the image in the live window. Do not forget to readjust the focus of the image on the screen. (Use Control + mouse scroll to adjust the focus)
- 3. Turn on Range Indicator at the bottom allows you to see whether there is overexposure. <u>Avoid overexposure!!!</u>
- 4. Change the exposure time to adjust the brightness of the image. Then click Stop to stop the live view
- 5. Then click on a different channel to view and adjust exposure time for that channel. Do this for all the channels you selected.







Acquire a multichannel image

- 1. After adjusting the exposure times for all the channels, click on Stop to stop the live view.
- 2. Then click Snap to capture the image with the channels you've selected.
- 3. After acquisiton, you will see a merged channel image. To view single channel instead of multiple channels simultaneously, check the single channel option, and click on the corresponding channel you wish to view.





Adjust image brightness and background post-acquisition

- 1. After image acquisition, you can adjust the image brightness through the histogram at the bottom of the software.
- 2. On the histogram, the right vertical bar moving to the left increases the intensity. The left vertical bar moving to the right darkens the background.





Acquire apotome images

1. Push the Apotome slider completely in to activate the apotome. Pull the slider out at the first stop or pull it out fully if you do not need to use the apotome.

- 2. You can acquire apotome images (Large scan) by checking the ApoTome High Grid or ApoTome Medium Grid acquisition settings.
- 3. For our 40x-100x objectives, use AptoTome High Grid settings, otherwise the outcome image is not good. (Very important!!)
- 4. Also make sure you pushed in the apotome slider.



👻 🏢 ApoTome Acqu	Jisition	
ApoTome Active		
Channel Calibrated		
Phase Images	5	•
Live Mode	Grid Visible	•



Acquire apotome images

5. After an apotome image is acquired, you can adjust the results through the Apotome options listed below the captured image. Phase errors are often used to get rid of grids.



-	Display ApoTome		
	Display mode	Optical sectioning	Normalization
it	Enable correction		
	Correction	Phase errors)
n	Fourier filter	Phase errors	Grid 428.99 L/mm
	Deconvolution	Local bleaching	
		Global bleaching	
5			
21		_	
	Create image		

Acquire a z-stack image

- 1. To enable Z-stack imaging, check the Z-stack box to activate the Z-stack window
- 2. Before defining the Z range, first adjust the exposure time for the brightest focal plane to make sure the channels do not overexpose at this focal plane.
- 3. To define the top and bottom of the z range, stay on live view, press Ctrl, and scroll the mouse wheel to change the focal plane. When you change the focal plane, you will find the real time z-position in the Z-stack window changes as you go.
- 4. First find the bottom z position of the cell/tissue and click Set First to record this position. Then scroll the mouse wheel the opposite direction and find the top z position of the cell/tissue, then click Set Last to record this position.
- 5. After the z range is defined, to change number of slices, first click Keep Slice, then type in the number of slices you prefer. To change Interval, click Keep Interval first, then type in the number for interval.
- 6. Then click Start Experiment to acquire the Z-stack image. Do not click Snap.







Acquire a Time-Series image

- 1. To enable time lapse imaging, check the Time Series box to activate the Time Series window. Uncheck others functions that you do not want to use. (Z-Stack, Tiles, Time Series)
- 2. In the Time Series window, set up the Duration (over all imaging time,) and the interval (time gap between each image.)
- 3. Then click Start Experiment.

Once you click Start Experiment, all the checked boxes will apply. So you can combine Z-stack with Time series and with Tiles.



" Multidimensional Acquisitio	n
👻 🕘 Time Series	🗸 Show All 🛛 🖉
Duration	10 Cycles V
As Long as Possible	
Interval	0.0 🗘 s 🔻
Use Camera Streaming if Possib	le
Use Burst Mode if Possible	
	1 Measure Speed
• Start	
Start Next Time Slice	
Stop	
Pause Begin	
Pause End	
Interval Time	
Set Marker	
Interactive Switches	
t Experiment Information	2

Acquire a Tile image

- 1. You can acquire tile images (Large scan) by checking the Tiles box.
- 2. Then you will find a new Tiles window and a Tiles section.
- 3. Then click Live to view the tile window simultaneously with the live window. (See next page)





- 4. Use the tool bar on top and the tool bar on the left Tiles functions to set up Large image scan area.
- 5. Use positions area to set up multi-points.
- 6. After an area is set, click Start Experiment to start acquisition.



Save and Export

- After image acquisition, save the image through File→Save. Save as .czi format.
- 2. DO NOT save anything on the desktop or C drive. Save everything in the Data drive (D drive) under your own folder.
- 3. Make sure you save your data as .cvi before exporting it to Tiff files.
- 4. Please copy your data when imaging done. Due to limited computer space, the core does not store image data for users.

ZEN	C:\User	rs\cshar	ng\Pictures\t	emp\Snap-1	104.czi ·	- ZEN 3
File	Edit	View	Acquisition	Graphics	Tools	Wind
	New Do	ocument	t	Ctrl+Shift	+N	
	New Im	age		Ctrl	+N	T
È	Open			Ctrl	+0	
	Save			Ctr	1+S	Арр
a	Save As	CZI		Ctrl+Shif	t+S	
	Save As	with Op	otions			
	Rename	2				(
	Delete					
	Export/I	Import			•	
	Send to	arivis V	ision4D			us
۲	Send to	ZEN bla	ack			
84.	Add to 2	ZEN Co	nnect Project			
	Close			Ctrl-	+F4	×
	Save All					

Save and Export

- 1. After saving the original image data. Export as TIF file through File→Export/Import→ Export.
- 2. This will bring you to the export window.
- 3. Select Tiff, Merged channels image, Individual channel images and use channel names.
- 4. Select the correct folder for exported images.
- 5. Then click Apply to export.

File	Edit	View	Acquisition	Graphics	Tools	Window	Help	
	New Do	ocument	t	Ctrl+Shift	+N			
	New Im	age		Ctrl	+N			()///4/
Đ	Open			Ctrl	+0	- <u>8</u> -2	•	Live.czi 🛛 Tile
周	Save			Ctr	l+S	Applicatio	ons	
믭	Save As	CZI		Ctrl+Shift	t+S	-	*	
	Save As	with O	ptions					
	Rename	:				45 F	Reuse	20
	Delete							
	Export/I	Import			•	Quick Expo	ort	
	Send to	arivis V	ision4D			Export		Ctrl+6
۲	Send to	ZEN bl	ack			Movie Exp	ort	
84.	Add to 2	ZEN Co	nnect Project			OME TIFF	Export	
×	Close			Ctrl+	+F4	ZVI Export		

Locate Acquisition Processing Analysis Applications
Function: Image Export
Single Batch Apply
Method
Parameters Show All
Settings 🔍 👻 🛪 🗸
Filetype Tiff Format (64Bit) (BigTiff)
■ BigTIFF ■ Tiff Tiles ■ Convert to 8 bit
Operation Lossless
Resize 1
Original data Shift Pixel
 Apply display curve and channel color Burn in graphics
Zoom 100 % 🗘 💿
 ✓ Merged channels image ✓ Individual channel images ✓ Use channel names
 Use full set of dimensions Define subset
Export To D:\Chaowei 🔚
✓ Create folder ☐ Generate xml file
Prefix Snap-104