

Year: _____

AED MONTHLY INSPECTION TRACKER

Item	Month (Initial in each box upon inspection)											
	Sept	Oct	Nov	Dec	Jan	Feb	Mar	April	May	June	July	Aug
General/Visual check: AED Clean												
All parts accounted for												
No parts past expiration												
Battery check												
Status indicator check												
Pads check- not expired												
Cabinet is secure												

Inspect the AED once every 30 days

This includes insuring that:

- The AED is located in its appropriate location
- The AED is clean and clear of any dirt or debris
- The AED has 1 set of sealed pads within the expiration date
- The AED's readiness indicator shows ready for use.