

How to use *neffy*



1 Remove

Remove *neffy* from packaging. Pull open the packaging to remove the *neffy* nasal spray device.



2 Hold

Hold device as shown. Hold the device with your thumb on the bottom of the plunger and a finger on either side of the nozzle.


- Do not pull or push on the plunger.
- Do not test or prime (pre-spray). Each device has only 1 spray.



3 Insert

Insert the nozzle into a nostril until your fingers touch your nose. Keep the nozzle straight into the nose pointed toward your forehead.

- Do not point (angle) the nozzle to the nasal septum (wall between your 2 nostrils) or outer wall of the nose.



4 Press

Press plunger up firmly until it snaps up and sprays liquid into the nostril.

- Do not sniff during or after the dose is given. If any liquid drips out of the nose, you may need to give a second dose of *neffy* after checking for symptoms.

Check for symptoms after the first dose of Neffy. If symptoms get worse or come back, give second dose of Neffy in the same nostril starting 5 minutes after the first dose

Get emergency medical help – call 9-1-1 – for further treatment of the allergic emergency (anaphylaxis) after using Neffy. Tell EMS personnel that Neffy nasal spray has been used and the time it was administered. Give EMS personnel the used Neffy device.